

Stroke-Fast-Track-API

Users

GET Get a user from username(admin)

http://localhost:3000/users/docter

HEADERS

User-Agent	insomnia/10.0.0
------------	-----------------

GET Me

http://localhost:3000/users/me

HEADERS

User-Agent	insomnia/10.0.0
------------	-----------------

GET Get Users

http://localhost:3000/users

HEADERS

User-Agent	insomnia/9.3.3
------------	----------------

POST Create a user

 This documentation is not published yet, this is only a preview. Changes to the publishing settings will reload this page.

HEADERS

Content-Type application/json

User-Agent insomnia/9.3.2

Body raw (json)

json

```
{
  "name": "Doctor Dome",
  "password" : "1234",
  "username" : "docter5",
  "email" : "docter5@hospital.com"
}
```

PATCH Update user

http://localhost:3000/users/update/docter4

HEADERS


Content-Type application/json

User-Agent insomnia/10.0.0

Body raw (json)

json

```
{
  "role": "Docter"
}
```

 This documentation is not published yet, this is only a preview. Changes to the publishing settings will reload this page.

localhost:3000/users/login

HEADERS

Content-Type application/json

User-Agent insomnia/9.3.2

Body raw (json)

json

```
{
  "usernameOrEmail": "admin",
  "password" : "1234"
}
```

POST Logout

http://localhost:3000/users/logout

HEADERS

User-Agent insomnia/9.3.2

PATCH Set Status (admin)

http://localhost:3000/users/setStatus

HEADERS

Content-Type application/json

 This documentation is not published yet, this is only a preview. Changes to the publishing settings will reload this page.

json

```
{  
  "username": "docter2",  
  "status": "true"  
}
```

DELETE Delete User

http://localhost:3000/users/delete

HEADERS

Content-Type application/json

User-Agent insomnia/9.3.3

Body raw (json)

json


```
{  
  "username": "docter5"  
}
```

Patient

GET Get Patient

http://localhost:3000/patient/

HEADERS

 This documentation is not published yet, this is only a preview. Changes to the publishing settings will reload this page.

POST Create Patient

http://localhost:3000/patient/create

HEADERS

Content-Type application/json

User-Agent insomnia/9.3.3

Body raw (json)

json

```
{
  "name": "John Doe",
  "gender": "Male",
  "dob": "1990-01-01",
  "weight": 75,
  "height": 180
}
```

PATCH Edit Patient

http://localhost:3000/patient/updte/9576f220-88d0-4677-9f6c-9b0d755d7054

HEADERS

Content-Type application/json

User-Agent insomnia/9.3.3

Body raw (json)

 This documentation is not published yet, this is only a preview. Changes to the publishing settings will reload this page.

```
{
  "name": "James",
  "gender": "Male",
  "dob": "2002-11-20T17:00:00.000Z",
  "weight": 75,
  "height": 183,
  "reg_id": "1234567890987"
}
```

DELETE Delete Patient

http://localhost:3000/patient/delete/921fb2c6-92a0-4161-96b7-01097d1665c9

HEADERS

Content-Type application/json

User-Agent insomnia/9.3.3

Body raw (json)

json

```
{
  "name": "John Doe",
  "gender": "Male",
  "age": 34,
  "dob": "1989-12-31T17:00:00.000Z",
  "weight": 75,
  "height": 180,
  "address": "",
  "symptoms": "ตาแดง"
}
```

Cases

GET Get case by patient id

GET Get a case by patient_id

 This documentation is not published yet, this is only a preview. Changes to the publishing settings will reload this page.

HEADERS

User-Agent	insomnia/10.0.0
------------	-----------------

GET Get all cases

http://localhost:3000/cases/

HEADERS

User-Agent	insomnia/9.3.3
------------	----------------

GET Get active cases

http://localhost:3000/cases/active/only

HEADERS

User-Agent	insomnia/10.0.0
------------	-----------------

GET Get admit cases

http://localhost:3000/cases/admit/only

HEADERS

User-Agent	insomnia/10.0.0
------------	-----------------

 This documentation is not published yet, this is only a preview. Changes to the publishing settings will reload this page.

HEADERS

Content-Type	application/json
User-Agent	insomnia/9.3.3

Body raw (json)

```
json

{
  "name": "James Doe",
  "gender": "Male",
  "dob": "1996/01/01",
  "weight": 80,
  "height": 178,
  "address": "456 Oak St",
  "onset": "2024-09-16T22:13:01.404+07:00"
}
```

PATCH Edit a case

http://localhost:3000/cases/update/9f5f5f58-e472-47de-9e6c-f7d0381654d6

HEADERS

Content-Type	application/json
User-Agent	insomnia/9.3.3

Body raw (json)

```
json

{
```



```
"status": "Admit"
```

 This documentation is not published yet, this is only a preview. Changes to the publishing settings will reload this page.

PATCH Finish a case

http://localhost:3000/cases/finish/9f5f5f58-e472-47de-9e6c-f7d0381654d6

HEADERS

User-Agent	insomnia/10.0.0
------------	-----------------

DELETE Delete a case

http://localhost:3000/cases/delete/cb0c133a-71e1-4483-9d0b-bcb35066520e

HEADERS

Content-Type	application/json
--------------	------------------

User-Agent	insomnia/9.3.3
------------	----------------

B.E.F.A.S.T

GET Get BEFAST by case_id

http://localhost:3000/cases/befast/bdc82dba-d865-4d96-ae9b-157d09006ba7

HEADERS

User-Agent	insomnia/10.0.0
------------	-----------------

POST Add BEFAST by case_id

POST Add BEFAST by case_id

 This documentation is not published yet, this is only a preview. Changes to the publishing settings will reload this page.

HEADERS

User-Agent	insomnia/10.0.0
------------	-----------------

PATCH Update BEFAST by case_id

http://localhost:3000/cases/update-befast/08310aaa-a800-4962-a67a-da6c67a03530

HEADERS

Content-Type	application/json
--------------	------------------

User-Agent	insomnia/10.0.0
------------	-----------------

Body raw (json)

json

```
{
  "b": "11",
  "b": "11",
  "e": "01",
  "f": "01",
  "a": "011",
  "s": "01",
  "t": "01"
}
```

CT Scan

GET Get CT Scan for a case

http://localhost:3000/cases/ct_result/123233c9-1f5d-4dba-b5d6-881acac7dade

HEADERS

 This documentation is not published yet, this is only a preview. Changes to the publishing settings will reload this page.

User-Agent	insomnia/10.0.0
------------	-----------------

POST Add CT Scan for a case

http://localhost:3000/cases/ct_result/01e599f5-8b47-443e-9e32-21e9ead1878c

HEADERS

User-Agent	insomnia/10.0.0
------------	-----------------

PATCH Edit CT Scan for a case

http://localhost:3000/cases/ct_result/01e599f5-8b47-443e-9e32-21e9ead1878c

HEADERS

Content-Type	application/json
--------------	------------------

User-Agent	insomnia/10.0.0
------------	-----------------

Body raw (json)

```
json
{
  "result": "Constricted"
}
```

NIHSS

GET Get NIHSS for a case

http://localhost:3000/cases/nihss/d603d610-c8b3-49cb-bef0-d0589fc5588d

 This documentation is not published yet, this is only a preview. Changes to the publishing settings will reload this page.

User-Agent	insomnia/10.0.0
------------	-----------------

POST Add NIHSS for a case

http://localhost:3000/cases/nihss/01e599f5-8b47-443e-9e32-21e9ead1878c

HEADERS

Content-Type	application/json
--------------	------------------

User-Agent	insomnia/10.0.0
------------	-----------------

Body raw (json)

```
json
```

```
{  
  "score": "10",  
  "round": "5"  
}
```

PATCH Edit NIHSS for a case

http://localhost:3000/cases/nihss/9f5f5f58-e472-47de-9e6c-f7d0381654d6

HEADERS

Content-Type	application/json
--------------	------------------

User-Agent	insomnia/10.0.0
------------	-----------------

Body raw (json)

```
{  
  "score": "12",  
  "round": "0"  
}
```

Thrombolytic

GET Get thrombolytic for a case

http://localhost:3000/cases/thrombolytic/d603d610-c8b3-49cb-bef0-d0589fc5588d

HEADERS

User-Agent	insomnia/10.0.0
------------	-----------------

POST Add for a case

http://localhost:3000/cases/thrombolytic/01e599f5-8b47-443e-9e32-21e9ead1878c

HEADERS

Content-Type	application/json
User-Agent	insomnia/10.0.0

PATCH Edit for a case

http://localhost:3000/cases/thrombolytic/01e599f5-8b47-443e-9e32-21e9ead1878c

HEADERS

 This documentation is not published yet, this is only a preview. Changes to the publishing settings will reload this page.

Injection

GET Get Injection

http://localhost:3000/cases/injection/d603d610-c8b3-49cb-bef0-d0589fc5588d

HEADERS

User-Agent	insomnia/10.0.0
------------	-----------------

POST Create Injection

http://localhost:3000/cases/injection/9f5f5f58-e472-47de-9e6c-f7d0381654d6

HEADERS

User-Agent	insomnia/10.0.0
------------	-----------------

PATCH Edit injection

http://localhost:3000/cases/injection/9f5f5f58-e472-47de-9e6c-f7d0381654d6

HEADERS

User-Agent	insomnia/10.0.0
------------	-----------------